

Anti-Npy4R Rabbit Polyclonal Antibody For Wb, Elisa - P50391

Item Number: CM0027128



Introduction

NPY4R Rabbit pAb for WB and ELISA, targeting human NPY4R. Recombinant immunogen (aa 1-80). Validated in human samples. 42kDa, shipped on ice. Ideal for neuropeptide receptor research.

[Learn More](#)

Parameter	Value
Product Name	NPY4R Rabbit pAb
Remarks/Alias	Y4; PP1; PPYR1; NPY4-R; NPY4R
Species	Human
GeneID	5540
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P50391
Protein Weight	42kDa
Shipping	Ice bag