

Anti-Npy5R Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q15761

Item Number: CM0011535

Introduction

High-quality rabbit monoclonal antibody against human NPY5R (SWISS: Q15761). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for neuropeptide Y receptor research.

[Learn More](#)

Parameter	Value
Product Name	NPY5R Rabbit mAb
Remarks/Alias	NPYR5; NPY5-R; NPY5-R; NPY5R
Species	Human
Gene ID (Human)	4889
Gene ID	4889
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NPY5R (Q15761).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15761
Protein Weight	51kDa
Shipping	Ice bag