

Anti-Nqo1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15559

Item Number: CM0014171



Introduction

Rabbit polyclonal antibody targeting human NQO1 (quinone reductase 1), validated for Western blot, IHC paraffin, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Suitable for investigating redox metabolism, chemoresistance, and proteasomal pathways.

[Learn More](#)

Parameter	Value
Product Name	NQO1 Rabbit pAb
Remarks/Alias	DTD; QR1; DHQU; DIA4; NMOR1; NMORI; NQO1
Species	Human
GeneID (Human)	1728
GeneID	1728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human NQO1 (NP_000894.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15559
Protein Weight	31kDa
Shipping	Ice bag