

Anti-Nqo1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P15559

Item Number: CM0006410



Introduction

Recombinant rabbit monoclonal NQO1 antibody for WB, IHC-P, IP, ELISA. Detects endogenous NQO1 (~31kDa) in human, mouse, and rat. Ideal for studying quinone metabolism, oxidative stress, and proteasomal regulation.

[Learn More](#)

Parameter	Value
Product Name	NQO1 Rabbit mAb
Remarks/Alias	DTD; QR1; DHQU; DIA4; NMOR1; NMORI; NQO1
Species	Human
GeneID (human)	1728
GeneID	1728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human NQO1 (NP_000894.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15559
Protein Weight	31kDa
Shipping	Ice bag