

Anti-Nqo1 Rabbit Monoclonal Antibody (Ko Validated) For Wb, If-P, Elisa - P15559

Item Number: CM0002891



Introduction

High-quality rabbit monoclonal antibody against NQO1, validated for knockout specificity. Suitable for western blot (WB), immunofluorescence (IF-P), and ELISA assays in human, mouse, and rat samples. Ideal for studying quinone detoxification, oxidative stress, and cancer-related pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NQO1 Rabbit mAb
Remarks/Alias	DTD; QR1; DHQU; DIA4; NMOR1; NMORI; [KD Validated] NQO1
Species	Human
Gene ID (Human)	1728
Gene ID	1728
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NQO1 (P15559).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15559
Protein Weight	31kDa
Shipping	Ice bag