

Anti-Nqo1 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, If-P, Elisa - P15559

Item Number: CM0007544



Introduction

KO-validated NQO1 rabbit monoclonal antibody for WB, IHC-P, IF-P, ELISA. Detects human, mouse, rat NQO1 (31 kDa). Ideal for redox biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NQO1 Rabbit mAb
Remarks/Alias	DTD; QR1; DHQU; DIA4; NMOR1; NMORI; [KD Validated] NQO1
Species	Human
Gene ID (Human)	1728
Gene ID	1728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human NQO1 (NP_000894.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15559
Protein Weight	31 kDa
Shipping	Ice bag