

## Anti-Nqo1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P15559

Item Number: CM0015624



### Introduction

Rabbit polyclonal antibody against human NQO1 (NAD(P)H dehydrogenase [quinone] 1). Suitable for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Used in redox biology and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | NQO1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | DTD; QR1; DHQU; DIA4; NMOR1; NMORI; NQO1                                                                        |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 1728                                                                                                            |
| GeneID           | 1728                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human NQO1 (NP_000894.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P15559                                                                                                          |
| Protein Weight   | 31kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |