

Anti-Nqo2 Monoclonal Antibody For Wb, Elisa - P16083

Item Number: CM0009511



Introduction

Rabbit monoclonal antibody against human NQO2 (Quinone reductase 2). Suitable for Western blot and ELISA; cross-reacts with mouse and rat. A useful tool for studying cellular detoxification pathways and vitamin K-dependent carboxylation.

[Learn More](#)

Parameter	Value
Product Name	NQO2 Rabbit mAb
Remarks/Alias	QR2; DHQV; DIA6; NMOR2; NQO2
Species	Human
Gene ID (Human)	4835
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NQO2 (P16083)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P16083
Protein Weight	26 kDa
Shipping	Ice bag