

Anti-Nqo2 Polyclonal Antibody For Wb, Elisa - P16083

Item Number: CM0014404



Introduction

NQO2 rabbit polyclonal antibody for WB and ELISA. Recognizes human and mouse NQO2 (QR2), a cytoplasmic quinone reductase involved in detoxification. Ideal for studies on vitamin K metabolism and detoxification pathways.

[Learn More](#)

Parameter	Value
Product Name	NQO2 Rabbit pAb
Remarks/Alias	QR2; DHQV; DIA6; NMOR2; NQO2
Species	Human
Gene ID (Human)	4835
Gene ID	4835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-231 of human NQO2 (NP_000895.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16083
Protein Weight	26kDa
Shipping	Ice bag