

Anti-Nqo2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P16083

Item Number: CM0016468



Introduction

Rabbit polyclonal antibody against human NQO2, validated for WB, IHC-P, ELISA. Suitable for quinone reductase studies in detoxification and vitamin K biosynthesis. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	NQO2 Rabbit pAb
Remarks/Alias	QR2; DHQV; DIA6; NMOR2; NQO2
Species	Human
GeneID (Human)	4835
GeneID	4835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-231 of human NQO2 (NP_000895.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16083
Protein Weight	26kDa
Shipping	Ice bag