


Anti-Nr1H4 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Ri1

Item Number: CM0018062



Introduction

High-quality rabbit polyclonal antibody against human NR1H4 (FXR), validated for Western blotting (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for bile acid receptor studies and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	NR1H4 Rabbit pAb
Remarks/Alias	BAR; FXR; HRR1; HRR-1; PFIC5; RIP14; NR1H4
Species	Human
Gene ID (Human)	9971
Gene ID	9971
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human NR1H4 (NP_001193922.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96Ri1
Protein Weight	56kDa
Shipping	Ice bag