


Anti-Nr1I3 Rabbit Polyclonal Antibody For Wb - Q14994

Item Number: CM0021277



Introduction

Rabbit polyclonal antibody against NR1I3 (CAR), suitable for WB, IF/ICC, IP, ELISA. Validated in human, mouse, rat. Ideal for studying nuclear receptor signaling and xenobiotic response.

[Learn More](#)

Parameter	Value
Product Name	NR1I3 Rabbit pAb
Remarks/Alias	CAR; CAR1; MB67; NR1I3
Species	Human
GeneID (Human)	9970
GeneID	9970
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 103-352 of human NR1I3 (NP_001070948.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14994
Protein Weight	40kDa
Shipping	Ice bag