

Anti-Nr2C1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13056

Item Number: CM0030206

Introduction

Rabbit polyclonal antibody against human NR2C1 (TR2), validated for WB, IF/ICC, ELISA; detects 67 kDa orphan nuclear receptor. Ideal for studying transcriptional repression, stem cells, and retinoic acid signaling.

[Learn More](#)

Parameter	Value
Product Name	NR2C1 Rabbit pAb
Remarks/Alias	TR2; NR2C1
Species	Human
Gene ID (Human)	7181
Gene ID	7181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human NR2C1 (NP_003288.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13056
Protein Weight	67kDa
Shipping	Ice bag