


Anti-Nr2C2 Rabbit Polyclonal Antibody For Wb, Elisa - P49116

Item Number: CM0030733



Introduction

Anti-NR2C2 rabbit polyclonal antibody validated for WB and ELISA; reacts with human, mouse, rat. Targets orphan nuclear receptor TR4, involved in transcription regulation, spermatogenesis, and embryonic development.

[Learn More](#)

Parameter	Value
Product Name	NR2C2 Rabbit pAb
Remarks/Alias	TR4; TAK1; NR2C2
Species	Human
GeneID (Human)	7182
GeneID	7182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human NR2C2 (NP_003289.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49116
Protein Weight	65kDa
Shipping	Ice bag