


Anti-Nr2C2/Tr4 Rabbit Monoclonal Antibody For Wb And Elisa - P49116

Item Number: CM0006395



Introduction

High-quality rabbit monoclonal antibody against human NR2C2/TR4, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying nuclear receptor signaling, spermatogenesis, and embryonic development.

[Learn More](#)

Parameter	Value
Product Name	NR2C2 / TR4 Rabbit mAb
Remarks/Alias	TR4; TAK1; NR2C2 / TR4
Species	Human
GeneID (Human)	7182
GeneID	7182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-555 of human NR2C2 / TR4 (NP_003289.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49116
Protein Weight	65kDa
Shipping	Ice bag