


Anti-Nr2E1/TLx Monoclonal Antibody For Wb And Elisa - Q9Y466

Item Number: CM0005742



Introduction

Rabbit monoclonal antibody against human NR2E1/TLX (nuclear receptor TLX). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Ideal for brain development and retinogenesis research.

[Learn More](#)

Parameter	Value
Product Name	NR2E1/TLX Rabbit mAb
Remarks/Alias	TLL; TLX; XTLL; NR2E1/TLX
Species	Human
GeneID (Human)	7101
GeneID	7101
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NR2E1 (Q9Y466).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y466
Protein Weight	43kDa
Shipping	Ice bag