


Anti-Nr2E1/Tlx Polyclonal Antibody For Wb, Elisa - Q9Y466

Item Number: CM0025228



Introduction

Rabbit polyclonal antibody against human NR2E1/TLX (Nuclear receptor TLX), validated for WB and ELISA, reacts with human and mouse. For studying brain development, retinal biology, and nuclear receptor signaling.

[Learn More](#)

Parameter	Specification
Product Name	NR2E1/TLX Rabbit pAb
Remarks/Alias	TLL; TLX; XTLL; NR2E1/TLX
Species	Human
GeneID (Human)	7101
GeneID	7101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 106-385 of human NR2E1/TLX (NP_003260.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y466
Protein Weight	43kDa
Shipping	Ice bag