

Anti-Nr2F1 Rabbit Polyclonal Antibody For IHC-P, Elisa - P10589

Item Number: CM0024592



Introduction

Rabbit polyclonal antibody against human NR2F1 (COUP-TF I), validated for IHC-P and ELISA. Cross-reactive with human. Useful for gene regulation and nuclear receptor studies. Molecular weight: 46kDa.

[Learn More](#)

Parameter	Value
Product Name	NR2F1 Rabbit pAb
Remarks/Alias	EAR3; BBOAS; EAR-3; SVP44; BBSOAS; ERBAL3; COUPTF1; TFCOUP1; COUP-TFI; TCFCOUP1; NR2F1
Species	Human
Gene ID (Human)	7025
Gene ID	7025
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 60-180 of human NR2F1 (NP_005645.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P10589
Protein Weight	46kDa
Shipping	Ice bag