


## Anti-Nr2F2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P24468

Item Number: CM0015251



### Introduction

Rabbit polyclonal antibody targeting human NR2F2 (COUP transcription factor 2). Suitable for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for research on transcription regulation and ovary development.

[Learn More](#)

| Parameter        | Value                                                                               |
|------------------|-------------------------------------------------------------------------------------|
| Product Name     | NR2F2 Rabbit pAb                                                                    |
| Remarks/Alias    | ARP1; ARP-1; CHTD4; NF-E3; SRXX5; SVP40; COUPTF2; COUPTFB; TFCOUP2; COUPTFII; NR2F2 |
| Species          | Human                                                                               |
| GeneID (Human)   | 7026                                                                                |
| GeneID           | 7026                                                                                |
| Immunogen        | Synthetic peptide corresponding to amino acids 1-100 of human NR2F2 (NP_066285.1)   |
| Source           | Rabbit                                                                              |
| Category         | Polyclonal Antibodies                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                   |
| SWISS            | P24468                                                                              |
| Protein Weight   | 46kDa                                                                               |
| Shipping         | Ice bag                                                                             |