


Anti-Nr5A2 Polyclonal Antibody For Wb, Ihc-P - O00482

Item Number: CM0018409



Introduction

Rabbit polyclonal NR5A2 antibody validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Essential for liver development, cholesterol metabolism, and hepatitis B virus research.

[Learn More](#)

Parameter	Value
Product Name	NR5A2 Rabbit pAb
Remarks/Alias	B1F, CPF, FTF, B1F2, LRH1, LRH-1, FTZ-F1, hB1F-2, FTZ-F1beta, NR5A2
Species	Human
GeneID (Human)	2494
GeneID	2494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-541 of human NR5A2 (NP_995582.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00482
Protein Weight	61kDa
Shipping	Ice bag