

Anti-Nras Rabbit Polyclonal Antibody For Wb - P01111

Item Number: CM0028193



Introduction

NRAS Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, and ELISA. Detects human, mouse, rat NRAS. Suitable for Ras-MAPK pathway and oncology studies. UniProt P01111.

[Learn More](#)

Parameter	Value
Product Name	NRAS Rabbit pAb
Remarks/Alias	NS6; CMNS; KRAS; NCMS; ALPS4; N-ras; NRAS1; NRAS
Species	Human
Gene ID (Human)	4893
Gene ID	4893
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 90-189 of human NRAS (NP_002515.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01111
Protein Weight	21kDa
Shipping	Ice bag