


Anti-Nrcam Polyclonal Antibody For Wb And Elisa - Q92823

Item Number: CM0022845



Introduction

Anti-NRCAM polyclonal antibody (rabbit) for WB and ELISA. Detects human and mouse NRCAM, a neuronal cell adhesion molecule critical for neurite outgrowth and node of Ranvier formation.

[Learn More](#)

Parameter	Value
Product Name	NRCAM Rabbit pAb
Remarks/Alias	NEDNMS; NRCAM
Species	Human
Gene ID (Human)	4897
Gene ID	4897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-115 of human NRCAM (Q92823).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92823
Protein Weight	144kDa
Shipping	Ice bag