

Anti-Nrd1 Rabbit Polyclonal Antibody For Wb, Elisa - O43847

Item Number: CM0029928

Introduction

Rabbit polyclonal antibody targeting human NRD1 (nardilysin), validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for studies on axonal maturation, myelination, and protease function.

[Learn More](#)

Parameter	Value
Product Name	NRD1 Rabbit pAb
Remarks/Alias	NRD1; hNRD1; hNRD2
Species	Human
Gene ID (Human)	4898
Gene ID	4898
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-400 of human NRD1 (NP_001095132.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43847
Protein Weight	132kDa
Shipping	Ice bag