

Anti-Nrf1 Rabbit Polyclonal Antibody For Wb, Elisa And Chip - Q16656

Item Number: CM0027256



Introduction

NRF1 Rabbit Polyclonal Antibody validated for WB, ELISA, and ChIP. Detects human, mouse, rat NRF1. Recombinant immunogen (aa 240-449). Suitable for transcription factor and mitochondrial biology research.

[Learn More](#)

Parameter	Value
Product Name	NRF1 Rabbit pAb
Remarks/Alias	ALPHA-PAL; NRF1
Species	Human
GeneID	4899
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-449 of human NRF1 (NP_001035199.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16656
Protein Weight	54kDa
Shipping	Ice bag