


Anti-Nrf2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16236

Item Number: CM0015730



Introduction

Rabbit polyclonal antibody to human NRF2 (NFE2L2), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for oxidative stress and cytoprotection research. SWISS: Q16236.

[Learn More](#)

Parameter	Value
Product Name	NRF2 Rabbit pAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH
Species	Human
GeneID (Human)	4780
GeneID	4780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 243-400 of human NRF2 (NP_006155.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16236
Protein Weight	68kDa
Shipping	Ice bag