


Anti-Nrf2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16236

Item Number: CM0014692



Introduction

High-quality rabbit polyclonal antibody against human NRF2, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for oxidative stress, cytoprotection, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	NRF2 Rabbit pAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH
Species	Human
GeneID (Human)	4780
GeneID	4780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 366-605 of human NRF2 (NP_006155.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16236
Protein Weight	68kDa
Shipping	Ice bag