


Anti-Nrf2 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q16236

Item Number: CM0014689



Introduction

Rabbit polyclonal antibody against human NRF2 (NP_006155.2, Q16236) for WB, IF-P, IHC-P, ELISA. Validated in human, mouse, rat. Suitable for oxidative stress, autophagy, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	NRF2 Rabbit pAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH
Species	Human
GeneID (Human)	4780
GeneID	4780
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16236
Protein Weight	68kDa
Shipping	Ice bag