


Anti-Nrg1 Polyclonal Antibody For Wb, Elisa - Q02297

Item Number: CM0017200



Introduction

NRG1 Rabbit pAb for Western blot and ELISA applications, reactive with human, mouse, and rat. Explore NRG1-ERBB signaling in cell growth, differentiation, and neurodevelopment. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	NRG1 Rabbit pAb
Remarks/Alias	GGF, HGL, HRG, NDF, ARIA, GGF2, HRG1, HRGA, SMDF, MST131, MSTP131, NRG1-IT2, NRG1
Species	Human
GeneID (Human)	3084
GeneID	3084
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human NRG1 (NP_039252.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02297
Protein Weight	70 kDa
Shipping	Ice bag