

Anti-Nrg4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wwg1

Item Number: CM0025628

Introduction

Rabbit polyclonal antibody to human neuregulin-4 (NRG4). Suitable for Western blot, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat. Research tool for ERBB tyrosine kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	NRG4 Rabbit pAb
Remarks/Alias	HRG4; NRG4
Species	Human
GeneID (Human)	145957
GeneID	145957
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human NRG4 (NP_612640.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WWG1
Protein Weight	13kDa
Shipping	Ice bag