


Anti-Nrip1 Polyclonal Antibody For Wb, Elisa - P48552

Item Number: CM0025175



Introduction

Rabbit polyclonal antibody targeting NRIP1 (RIP140), suitable for Western blot and ELISA applications. Detects human and mouse NRIP1, a transcriptional coregulator involved in circadian rhythm and renal development.

[Learn More](#)

Parameter	Value
Product Name	NRIP1 Rabbit pAb
Remarks/Alias	CAKUT3; RIP140; NRIP1
Species	Human
GeneID (human)	8204
GeneID	8204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 939-1158 of human NRIP1 (NP_003480.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48552
Protein Weight	127kDa
Shipping	Ice bag