

Anti-Nrip2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bqi9

Item Number: CM0016869

Introduction

High-quality rabbit polyclonal antibody against human NRIP2 (Nuclear receptor-interacting protein 2), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	NRIP2 Rabbit pAb
Remarks/Alias	NRIP2
Species	Human
Gene ID (Human)	83714
Gene ID	83714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-281 of human NRIP2 (NP_113662.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQI9
Protein Weight	31kDa
Shipping	Ice bag