


Anti-Nrxn1 Rabbit Polyclonal Antibody For Wb, Elisa - P58400 / Q9ULb1

Item Number: CM0012933



Introduction

High-quality anti-NRXN1 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse neurexin-1-beta, a key synaptic adhesion molecule involved in synapse formation. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	NRXN1 Rabbit pAb
Remarks/Alias	PTHSL2; SCZD17; Hs.22998; NRXN1
Species	Human
GeneID (Human)	9378
GeneID	9378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-310 of human NRXN1 (NP_001129131.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P58400 / Q9ULB1
Protein Weight	15kDa / 46kDa / 161kDa / 164kDa / 169kDa
Shipping	Ice bag