

Anti-Nrxn3 Polyclonal Antibody For Wb, Elisa - Q9Hdb5/Q9Y4C0

Item Number: CM0016018

Introduction

Rabbit polyclonal antibody against human NRXN3, reactive with human, mouse, and rat. Validated for WB and ELISA. Ideal for studying neurexin-3's role in synaptic adhesion and receptor modulation.

[Learn More](#)

Parameter	Value
Product Name	NRXN3 Rabbit pAb
Remarks/Alias	C14orf60; NRXN3
Species	Human
GeneID (human)	9369
GeneID	9369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 373-432 of human NRXN3 (NP_620426.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HDB5/Q9Y4C0
Protein Weight	43kDa/47kDa/50kDa/69kDa/117kDa/153kDa/180kDa
Shipping	Ice bag