


Anti-Nsa2 Rabbit Polyclonal Antibody For Wb And Elisa - O95478

Item Number: CM0015447



Introduction

High-quality rabbit polyclonal antibody against human NSA2 (TINP1/HCLG1), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	NSA2 Rabbit pAb
Remarks/Alias	HCLG1; TINP1; CDK105; HCL-G1; HUSSY29; HUSSY-29; NSA2
Species	Human
Gene ID (Human)	10412
Gene ID	10412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human NSA2 (NP_055701.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95478
Protein Weight	30kDa
Shipping	Ice bag