


Anti-Nsd3/Whsc1L1 Polyclonal Antibody For Wb, Elisa - Q9Bz95

Item Number: CM0019754



Introduction

Rabbit polyclonal antibody targeting human NSD3/WHSC1L1, validated for WB and ELISA. Reacts with human and mouse. Useful for studying H3K4me2 and H3K27me2 epigenetic marks in chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	NSD3/WHSC1L1 Rabbit pAb
Remarks/Alias	KMT3F; KMT3G; WHISTLE; WHSC1L1; pp14328; NSD3/WHSC1L1
Species	Human
Gene ID	54904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-590 of human NSD3/WHSC1L1 (NP_075447.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZ95
Protein Weight	162kDa
Shipping	Ice bag