

Anti-Nsf Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P46459

Item Number: CM0016897



Introduction

High-quality rabbit polyclonal antibody targeting human NSF (vesicle-fusing ATPase) validated for WB, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat samples. Suitable for vesicle trafficking and Golgi function studies.

[Learn More](#)

Parameter	Value
Product Name	NSF Rabbit pAb
Remarks/Alias	SKD2; DEE96; SEC18; NSF
Species	Human
Gene ID (Human)	4905
Gene ID	4905
Immunogen	Recombinant fusion protein containing amino acids 251-400 of human NSF (NP_006169.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46459
Protein Weight	83kDa
Shipping	Ice bag