

Anti-Nsfl1C Rabbit Polyclonal Antibody (P47) For Wb, If/Icc, Ip, Elisa - Q9Unz2

Item Number: CM0025325



Introduction

High-quality rabbit polyclonal antibody against human NSFL1C/p47, validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for Golgi, mitosis, and VCP-related research.

[Learn More](#)

Parameter	Value
Product Name	NSFL1C Rabbit pAb
Remarks/Alias	P47; UBX1; UBXD10; UBXN2C; dj776F14.1
Species	Human
Gene ID (Human)	55968
Gene ID	55968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human NSFL1C (NP_057227.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNZ2
Protein Weight	41kDa
Shipping	Ice bag