

Anti-NSL1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Iy1

Item Number: CM0025783



Introduction

NSL1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human kinetochore-associated protein NSL1, crucial for mitotic chromosome alignment and segregation. Ideal for cell cycle and kinetochore research.

[Learn More](#)

Parameter	Value
Product Name	NSL1 Rabbit pAb
Remarks/Alias	DC8, MIS14, C1orf48, NSL1
Species	Human
GeneID (human)	25936
GeneID	25936
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human NSL1 (NP_056286.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96IY1
Protein Weight	32kDa
Shipping	Ice bag