

Anti-Nsmf Polyclonal Antibody For Wb, Elisa - Q6X4W1

Item Number: CM0027837

Introduction

NSMF rabbit polyclonal antibody for WB and ELISA; detects human and mouse NSMF (also known as NELF or HH9), a 60 kDa NMDA receptor synaptonuclear signaling factor involved in neuronal migration and dendrite cytoarchitecture. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	NSMF Rabbit pAb
Remarks/Alias	HH9; NELF; NSMF
Species	Human
GeneID (Human)	26012
GeneID	26012
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 226-505 of human NSMF (NP_001124443.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6X4W1
Protein Weight	60kDa
Shipping	Ice bag