

Anti-Nsrp1 Polyclonal Antibody For Wb, Elisa - Q9H0G5

Item Number: CM0022667

Introduction

Rabbit polyclonal antibody targeting NSRP1 (CCDC55/NSrp70), validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for studies on pre-mRNA splicing and nuclear speckle dynamics.

[Learn More](#)

Parameter	Specification
Product Name	NSRP1 Rabbit pAb
Remarks/Alias	CCDC55; NSrp70; HSPC095; NEDSSBA; NSRP1
Species	Human
GeneID (Human)	84081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NSRP1 (NP_115517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0G5
Protein Weight	66kDa
Shipping	Ice bag