

Anti-Nsun2 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q08J23

Item Number: CM0006956

Introduction

NSUN2 Rabbit monoclonal antibody validated for Western blot, immunohistochemistry, immunofluorescence, immunoprecipitation, and ELISA. Detects human and mouse NSUN2, an RNA methyltransferase involved in tRNA/mRNA modification and development.

[Learn More](#)

Parameter	Value
Product Name	NSUN2 Rabbit mAb
Remarks/Alias	MISU; MRT5; SAKI; TRM4; NSUN2
Species	Human
Gene ID (Human)	54888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 617-708 of human NSUN2 (NP_060225.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08J23
Protein Weight	86 kDa
Shipping	Ice bag