

Anti-Nsun3 Rabbit Polyclonal Antibody For Wb & Elisa - Q9H649

Item Number: CM0016153

Introduction

Rabbit polyclonal antibody against human NSUN3, a mitochondrial tRNA methyltransferase associated with COXPD48. Validated for WB and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	NSUN3 Rabbit pAb
Remarks/Alias	MST077; COXPD48; MSTP077; NSUN3
Species	Human
Gene ID (Human)	63899
Gene ID	63899
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human NSUN3 (NP_071355.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H649
Protein Weight	38kDa
Shipping	Ice bag