


Anti-Nsun5 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96P11

Item Number: CM0025617



Introduction

Rabbit polyclonal antibody targeting human NSUN5 (NOL1/p120), validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Investigate the role of 28S rRNA methyltransferase in protein translation, development, and Williams-Beuren syndrome.

[Learn More](#)

Parameter	Value
Product Name	NSUN5 Rabbit pAb
Remarks/Alias	NOL1; p120; NOL1R; NSUN5A; WBSCR20; WBSCR20A; p120(NOL1); NSUN5
Species	Human
GeneID (Human)	55695
GeneID	55695
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-429 of human NSUN5 (NP_060514.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P11
Protein Weight	47kDa
Shipping	Ice bag