

Anti-Nt5E/Cd73 Rabbit Monoclonal Antibody - P21589

Item Number: CM0006191



Introduction

Recombinant rabbit monoclonal antibody against human NT5E/CD73, validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for nucleotidase activity and cancer research.

[Learn More](#)

Parameter	Value
Product Name	NT5E/CD73 Rabbit mAb
Remarks/Alias	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA
Species	Human
GeneID (Human)	4907
GeneID	4907
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 27-547 of human NT5E/CD73 (NP_002517.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21589
Protein Weight	57kDa/63kDa
Shipping	Ice bag