


Anti-Ntmt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bv86

Item Number: CM0029462



Introduction

High-quality rabbit polyclonal antibody against human NTMT1 (N-terminal methyltransferase). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for chromatin and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	NTMT1 Rabbit pAb
Remarks/Alias	NRMT; NRMT1; NTM1A; AD-003; HONT1A; C9orf32; METTL11A; NTMT1
Species	Human
GeneID (human)	28989
GeneID	28989
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV86
Protein Weight	25kDa
Shipping	Ice bag