


Anti-Ntpcr Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bsd7

Item Number: CM0018617



Introduction

NTPCR rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse NTPCR, a cancer-related nucleoside-triphosphatase. Ideal for nucleotide metabolism and oncology research.

[Learn More](#)

Parameter	Value
Product Name	NTPCR Rabbit pAb
Remarks/Alias	THEP1, C1orf57, HCR-NTPase, NTPCR
Species	Human
GeneID (Human)	84284
GeneID	84284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 91-190 of human NTPCR (NP_115700.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BSD7
Protein Weight	21kDa
Shipping	Ice bag