

Anti-Ntrk1/Ntrk2/Ntrk3 Polyclonal Antibody For Wb, If/Icc, Elisa - P04629/Q16620/Q16288

Item Number: CM0020772



Introduction

High-quality rabbit polyclonal antibody targeting NTRK1, NTRK2, NTRK3. Validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for neurotrophin signaling research and neurobiology assay development.

[Learn More](#)

Parameter	Specification
Product Name	NTRK1/NTRK2/NTRK3 Rabbit pAb
Species	Human
GeneID (Human)	4914/4915/4916
GeneID	4914/4915/4916
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-796 of human NTRK1/NTRK2/NTRK3 (NP_002520.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04629/Q16620/Q16288
Shipping	Ice bag