

Anti-Nts Polyclonal Antibody For Wb, If-P, Elisa - P30990

Item Number: CM0015218



Introduction

High-quality rabbit polyclonal antibody against human NTS, validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for studying neurotensin function in fat metabolism and smooth muscle contraction.

[Learn More](#)

Parameter	Value
Product Name	NTS Rabbit pAb
Remarks/Alias	NN; NT; NT/N; NTS1; NMN-125; NTS
Species	Human
GeneID (Human)	4922
GeneID	4922
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-148 of human NTS (NP_006174.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P30990
Protein Weight	20kDa
Shipping	Ice bag