

Anti-Nts Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P30990

Item Number: CM0013752



Introduction

NTS Rabbit pAb is a polyclonal antibody recognizing human, mouse, and rat neurotensin/neuromedin N. Validated for WB, IHC-P, and ELISA, it is ideal for exploring neurotensin's roles in fat metabolism and smooth muscle regulation.

[Learn More](#)

Parameter	Value
Product Name	NTS Rabbit pAb
Remarks/Alias	NN; NT; NT/N; NTS1; NMN-125; NTS
Species	Human
Gene ID (Human)	4922
Gene ID	4922
Immunogen	synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-170 of human NTS (NP_006174.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30990
Protein Weight	20kDa
Shipping	Ice bag