

Anti-Ntsr1 Polyclonal Antibody For Wb And Elisa - P30989

Item Number: CM0028068



Introduction

High-affinity neurotensin receptor type 1 (NTSR1) rabbit polyclonal antibody validated for WB and ELISA, reactive with human and mouse. Suitable for GPCR signaling and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	NTSR1 Rabbit pAb
Remarks/Alias	NTR; NTSR1
Species	Human
Gene ID (Human)	4923
Gene ID	4923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human NTSR1 (NP_002522.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30989
Protein Weight	46kDa
Shipping	Ice bag